

# MACHEREY-NAGEL

## CHROMAFIL<sup>®</sup> Xtra

### Syringe filters

**UPDATE 2019**  
**DUE TO THE BIG SUCCESS**  
**OF LAST YEAR'S CAMPAIGN**

Ask for  
**SPECIAL**  
**DISCOUNT!**  
**%**



For reliable and successful filtrations

- PET syringe filters
- PTFE syringe filters
- PES syringe filters
- Polyamide / Nylon syringe filters



**MACHEREY-NAGEL**

[www.mn-net.com](http://www.mn-net.com)



Chromatography

# CHROMAFIL® Xtra

## CHROMAFIL® Xtra syringe filters

CHROMAFIL® syringe filters are used for the filtration of suspended matter from liquid samples (1–100 mL) or gases. With CHROMAFIL®, rapid purification and removal of particles is very simple: just place the filter on the syringe, and you are ready for filtration. Contamination of sensitive instrumentation by solid impurities can be avoided, thus increasing lifetime of chromatographic columns and equipment. The filters can be used for sample preparation in HPLC, GC, ICP, AAS, TOC, DOC, IR, NMR, photometry, spectroscopy, etc.



Product	Membrane pore size	Membrane dia.	Pack of	REF
<b>Polyester (PET), for polar as well as nonpolar solvents</b>				
CHROMAFIL® Xtra PET-20/13	0.20 µm	13 mm	100	729222
CHROMAFIL® Xtra PET-20/25	0.20 µm	25 mm	100	729221
CHROMAFIL® Xtra PET-45/13	0.45 µm	13 mm	100	729223
CHROMAFIL® Xtra PET-45/25	0.45 µm	25 mm	100	729220
<b>Polytetrafluoroethylene (PTFE), for nonpolar liquids and gases</b>				
CHROMAFIL® Xtra PTFE-20/13	0.20 µm	13 mm	100	729208
CHROMAFIL® Xtra PTFE-20/25	0.20 µm	25 mm	100	729207
CHROMAFIL® Xtra PTFE-45/13	0.45 µm	13 mm	100	729209
CHROMAFIL® Xtra PTFE-45/25	0.45 µm	25 mm	100	729205
<b>Polyethersulfone (PES), for aqueous and slightly organic liquids</b>				
CHROMAFIL® Xtra PES-20/25	0.20 µm	25 mm	100	729240
CHROMAFIL® Xtra PES-45/25	0.45 µm	25 mm	100	729241
<b>Polyamide / Nylon (PA), for aqueous organic and aqueous medium polar liquids</b>				
CHROMAFIL® Xtra PA-20/13	0.20 µm	13 mm	100	729248
CHROMAFIL® Xtra PA-20/25	0.20 µm	25 mm	100	729212
CHROMAFIL® Xtra PA-45/13	0.45 µm	13 mm	100	729249
CHROMAFIL® Xtra PA-45/25	0.45 µm	25 mm	100	729213
<b>Special filters for ion chromatography (IC), for aqueous liquids</b>				
CHROMAFIL® Xtra IC-45/25	0.45 µm	25 mm	100	729258



Str. Dobrogeanu Gherea nr. 39  
Targu Mures, jud. Mures, Romania  
tel. 0365-421786  
tel. 0745110534  
fax. 0365-814199  
office@bioaqua.ro



www.mn-net.com

**MACHEREY-NAGEL**



MACHEREY-NAGEL GmbH & Co. KG · Neumann-Neander-Str. 6–8 · 52355 Düren · Germany

DE / International:

Tel.: +49 24 21 969-0

Fax: +49 24 21 969-199

E-mail: info@mn-net.com

CH:

Tel.: +41 62 388 55 00

Fax: +41 62 388 55 05

E-mail: sales-ch@mn-net.com

FR:

Tel.: +33 388 68 22 68

Fax: +33 388 51 76 88

E-mail: sales-fr@mn-net.com

US:

Tel.: +1 484 821 0984

Fax: +1 484 821 1272

E-mail: sales-us@mn-net.com

